

**INTENDED USE:**

**NATtrol™ Respiratory Panel Control** (NATRSPC-DIA) is a qualitative external control intended to be used with molecular assays that detect nucleic acids from the organisms listed in Table 1. The control can be used to monitor the steps of extraction, amplification, and detection, and also for training, and evaluation of operator variation.

**PRODUCT DESCRIPTION:**

NATtrol™ Respiratory Panel Control contains purified, intact bacteria and viruses that have been chemically modified to render them non-infectious and refrigerator stable. The control is supplied in a proprietary matrix.

**MATERIALS SUPPLIED:**

- 3 x 0.35 mL vials of Respiratory Panel Control 1 and 3 x 0.35 mL vials of Respiratory Panel Control 2.

**WARNINGS AND PRECAUTIONS:**

- NATtrol™ inactivation was carried out on virus and microorganism stocks used to formulate the control. The inactivation was verified in a standard microbiological growth protocol.
- This product contains inactivated viruses and microorganisms and materials of human and animal origin. Safe practices suggest that the controls be considered potentially infectious and to use Universal Precautions when handling.
- Refer to CDC guidelines and local regulations for handling and disposal.
- The matrix used in the manufacture of this product is treated with 0.09% sodium azide. It was manufactured from Human Serum Albumin that has been tested and found to be non-reactive at the donor level for HIV-1/HIV-2 Antibody, HBsAg and HCV Antibody by FDA licensed donor screening test methods. All materials are also tested for HIV-1 and HCV by FDA approved Nucleic Acid Test (NAT) methods.
- Heat inactivated Fetal Bovine Serum used in the manufacture of this product meets applicable USDA requirements for abattoir sourced animals, traceability, and country of origin. The materials were collected at USDA licensed establishments or legally imported from countries recognized by the USDA as negligible or controlled for risk for Bovine Spongiform Encephalopathy (BSE) and other exotic disease agents. Donor animals were inspected ante and postmortem at the abattoir as required by the USDA.
- Do not use past the expiration date on the label.
- To avoid cross-contamination, use separate pipette tips for all materials.

**RECOMMENDED STORAGE:**

- NATtrol™ Respiratory Panel Control should be stored at 2-8°C.

**INSTRUCTIONS FOR USE:**

- Mix tube vigorously for at least 5 secs.
- Process according to manufacturer's instructions for sample to result assays.
- Extract nucleic acid prior to use in downstream assays that are not sample to result.
- Each vial is intended for single use.

**LIMITATIONS:**

- **FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES**
- Quality control materials should be used in accordance with local, state, federal, and accreditation requirements.
- This product is not intended to replace the manufacturer's controls provided with the assay.

**EXPECTED RESULTS:**

The table below is for informational purposes, each laboratory must evaluate the control and establish their own performance criteria. NATrol™ Respiratory Panel Control is a qualitative control and does not have assigned analyte values or ranges.








**TABLE 1: Organisms/Viruses in Respiratory Panel Control**

Respiratory Panel Control	Organism/Virus (Strain)
<b>Respiratory Panel Control 1</b>	Adenovirus Type 3
	<i>B. parapertussis</i> (A747)
	<i>C. pneumoniae</i> (IOL-207)
	Coronavirus (OC43)
	Influenza B (B/Florida/02/06)
	Metapneumovirus 8 (Peru6-2003 <sup>1</sup> )
	Parainfluenza Type 1
	Parainfluenza Type 2
	RSV B (CH93(18)-18)
	SARS-Cov-2 - USA-WA1/2020 <sup>2</sup>
<b>Respiratory Panel Control 2</b>	<i>B. holmesii</i> (F061)
	<i>B. pertussis</i> (A639)
	Influenza A H1N1 (A/Singapore/63/04)
	Influenza A H3N2 (A/Brisbane/10/07)
	<i>M. pneumoniae</i> (M129)
	Parainfluenza Type 3
	Parainfluenza Type 4
	Rhinovirus 1A

<sup>1</sup>This product is sold by ZeptoMetrix under license from Vironovative B. V. under patent applications, including U.S. Patent Applications 10/371,099 and 10/371,12 and any patents that issue from applications related to PCT/NL02/00040 and PCT/US03/05271.

<sup>2</sup>This reagent was deposited by the Centers for Disease Control and Prevention and obtained through BEI Resources, NIAID, NIH: SARS-Related Coronavirus 2, Isolate USA-WA1/2020, NR-52281.

**Symbols used in the labeling of this product:**

	<b>Manufacturer</b>		<b>Temperature Limitation</b>
	<b>For Research Use Only</b>		<b>Expiration Date</b>
	<b>Catalog Number</b>		<b>Biological Risk</b>
	<b>Batch Code</b>		