

# Baby System<sup>®</sup>

## Graf's cradle (Rapid-Eco and Eco-Support)

**Ultrasound Positioning Device with arm guide for ultrasound probes - Medical speciality : medical sonographer/doctors (orthopedists, pediatricians)/nurses**

**Rapid-Eco** is the ultrasound positioning device which is designed to hold and position a newborn baby quickly and efficiently for a scan. On this piece of equipment the newborn is kept firmly still and in the correct position without needing the assistance of a third person, this is because the pressure from both sides of the equipment blocks the movement of the baby with a pressure which is in proportion to the baby's weight and hence the doctor is able to carry on the medical examination in a quick and precise way. **Rapid-Eco** is washable and made from materials of advanced technology. The Rapid-Eco gives doctors an easier and quicker way of working.



**Eco-Support** is the arm guide for ultrasound probes is made of steel and is fixed to the base of the Rapid-Eco. During the scan, the doctor is able to observe and analyses the data on the monitor, as the Eco-Support allows precise movement and hence avoids errors in the diagnosis.

Both products (**Rapid-Eco product code 5000** and **Eco-Support product code 5001**) can be purchased singularly upon request.

**DESIGNED, CERTIFIED AND AUTHORIZED BY Dr. Prof. Reinhard Graf**



### Additional Information

- Product code 5003 - HS 94029000 - size 27x51x80 cms - weight 7.3 KGS
- Alternative product : *Portable Graf's cradle product code 5002*

[CHECK PRICE ON LINE](#)